

COMeN Share with the World



NV6 Neonatal Ventilator

NV6 provides bubble CPAP ventilation. Bubbles vibrating in high frequency generated by the pressure generating device, generate high-frequency oscillations of the gas flow entering through the infant's respiratory system, which rapidly tap the child's lungs to improve respiratory function, reduce dead corners and carbon dioxide retention, promote sputum drainage and accelerate recovery.

NV6 also features high-flow nasal cannula (HFNC), which can achieve the performance as well as CPAP, and it is much more convenient to operate, with less damage to infant's nose.



Shenzhen Comen Medical Instruments Co., Ltd.

P/N: EN-NV6-4P-20200219-V1.0

Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Gongming Sub-district,Guangming New District, Shenzhen, 518106, P.R.China
Tel:+86-755-2640 8879 400-700-9488 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com

COMeN

Medical air oxygen mixer

- Stable input: At the air inlet end AIR / O₂ we adopted GENTEC gas connector and hose; the connector is made of high-quality copper alloy, and the surface is electroplated. A precise filtering device is installed internally to allow the clean gas source delivering steadily.
- Pressure imbalance alarm: In order to stabilize the pressure, KL-20 applies the advanced second order balance design and the principle of mechanical diaphragm balance. If the difference between oxygen and air pressure is too large, the differential pressure alarm will be triggered.
- Air-oxygen mixing: KL-20 adopts a high-precision proportioning valve. By adjusting the control knob, the accuracy of oxygen concentration output after air-oxygen mixing is controlled at an amazing 3%.
- Precise oxygen output: The GENTEC high-precision flowmeter is adopted at the oxygen output end of KL-20, which features aluminum alloy material, anodized surface and rear cut-off design, ensuring stable flow without being affected by outlet pressure.



High-quality CPAP generator

- Adopts imported nBubble PAP Valve.
- End-expiratory pressure setting range: 0-10cmH₂O.
- Automatic water level adjustment mechanism: maintain constant CPAP pressure level.
- Detachable water tank: drain excess condensate without interrupting CPAP.



Quiet medical air compressor (optional)

- Hierarchical noise cancellation technology: specifically for high and low frequency noise cancellation.
- Filtration technology: minimal 0.01um, adopts SMC polymer membrane dryer to filter out water molecules.



Servo-controlled humidifier

- Advanced humidification treatment system can restore the body's natural balance ,help build the body's physiological defense mechanism and make patients more comfortable.

