### COMON

## **COMEN** Share with the World



**C**€<sub>1639</sub>



Add: Floor 10, Floor 11 and Section C of Floor 12 of Building 1A &Floor 1 to Floor 5 of Building 2, FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district, Guangming District, Shenzhen, Guangdong, 518106,P.R. China Tel:+86-755-2640 8879 Fax:+86-755-2643 1232 Website: en.comen.com E-mail: info@szcomen.com

P/N: EN-S5-8P-20210528-V1.1

# S5 Defibrillator Monitor



#### Defibrillation

Manual defibrillation modes include synchronous cardioversion and asynchronous defibrillation.

#### Pacing

Having on-demand pacing and fixed pacing mode, for patients with cardiac arrest and acute severe slow arrhythmia, in vitro non-invasive time-saving and improve recovery success rate. emergency personnel in providing defibrillation and basic life support.

#### Monitor

5-lead ECG monitoring as standard, optional monitoring functions include SpO2, NIBP, PR and EtCO2 are also available for continuous monitoring of patient vital signs.

#### 3s

#### 3 steps

COMeN 01





#### **EN 1789 certificated**

S5 designed for ambulance and emergency transportation, after rigorous testing, it acquired EN 1789 certification.



#### **Compact yet lighter**

The S5 is perfect design of small size and lightweight(down to 5.6kg) defibrillator/monitor so as to be readily portable in the pre-hospital and hospital environment effectively.



#### **Longest lifespan**

Two kinds of battery capacity, both of them provide excellent performance, support fast charging up to 210 times of maximum discharges, satisfy clinical requirements.



#### **Better protection**

With the protection of IP44, S5 offers reliable performance in a variety of adverse conditions.

COMEN 02



# Advanced BTE Technology

Higher efficiency with Advanced Biphasic Truncated Exponential (BTE) waveform technology and automatic impedance compensation.

Up to 360J energy selection, for patients with high defibrillation thresholds such as myocardial infarction, obesity, high impedance, higher energy selection indicates greater defibrillation success rate.



# 1 for 2

Electrode paddles can be divided into large and small electrode paddles which can defibrillate for both adult and child respectively.

1 pair

Defibrillating operation can be completed on a pair of electrode paddles.